


Expersus Femoral Stem

BACKUP INFORMATION






Overview




Expersus Femoral Stem STD, 125°, cementless (ISO 5832-3 Ti6Al4V) with Head Ø28

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø28/Neck component selected				Stem Offset [mm] when Head Ø28 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1409	2	139,8	118,1	34,6	38,1	41,6	45,1	41,5	43,5	45,5	47,5
367-1410	3	145,0	123,2	35,5	39,0	42,5	46,0	43,5	45,5	47,5	50,0
367-1411	4	150,1	128,6	34,7	38,2	41,7	45,2	43,9	45,9	47,9	50,4
367-1412	5	154,1	132,6	34,9	38,4	41,9	45,4	44,7	46,7	48,7	51,2
367-1413	6	159,9	138,4	34,9	38,4	41,9	45,4	45,1	47,1	49,1	51,6
367-1414	7	164,9	143,5	35,0	38,5	42,0	45,5	45,9	47,9	49,9	52,4
367-1415	8	169,1	147,7	35,0	38,5	42,0	45,5	47,1	49,1	51,1	53,6
367-1416	9	180,0	158,6	35,2	38,7	42,2	45,7	48,4	50,4	52,4	54,9
367-1417	10	189,2	167,9	35,4	38,9	42,4	45,9	49,6	51,6	53,6	56,1




Expersus Femoral Stem STD, 125°, cementless (ISO 5832-3 Ti6Al4V) with Head Ø32

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø32/Neck component selected				Stem Offset [mm] when Head Ø32 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1409	2	139,8	118,1	34,1	38,1	42,1	46,1	41,2	43,5	45,8	48,1
367-1410	3	145,0	123,2	35,0	39,0	43,0	47,0	43,2	45,5	47,8	50,1
367-1411	4	150,1	128,6	34,2	38,2	42,2	46,2	43,6	45,9	48,2	50,5
367-1412	5	154,1	132,6	34,4	38,4	42,4	46,4	44,4	46,7	49,0	51,3
367-1413	6	159,9	138,4	34,4	38,4	42,4	46,4	44,8	47,1	49,4	51,7
367-1414	7	164,9	143,5	34,5	38,5	42,5	46,5	45,6	47,9	50,2	52,5
367-1415	8	169,1	147,7	34,5	38,5	42,5	46,5	46,8	49,1	51,4	53,7
367-1416	9	180,0	158,6	34,7	38,7	42,7	46,7	48,1	50,4	52,7	55,0
367-1417	10	189,2	167,9	34,9	38,9	42,9	46,9	49,3	51,6	53,9	56,2




Expersus Femoral Stem STD, 135°, cementless (ISO 5832-3 Ti6Al4V) with Head Ø28

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø28/Neck component selected				Stem Offset [mm] when Head Ø28 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1326	2	139,8	118,1	35,1	38,6	42,1	45,6	36,3	38,8	41,3	43,7
367-1327	3	145,0	123,2	35,6	39,1	42,6	46,1	37,5	40,0	42,5	44,9
367-1328	4	150,1	128,6	35,1	38,6	42,1	45,6	38,3	40,8	43,3	45,7
367-1329	5	154,0	132,6	35,1	38,6	42,1	45,6	38,9	41,4	43,9	46,3
367-1330	6	159,9	138,4	35,1	38,6	42,1	45,6	39,3	41,8	44,3	46,7
367-1331	7	164,9	143,5	35,1	38,6	42,1	45,6	40,0	42,5	45,0	47,4
367-1332	8	170,1	148,4	35,1	38,6	42,1	45,6	41,0	43,5	46,0	48,4
367-1333	9	180,0	158,6	35,1	38,6	42,1	45,6	42,1	44,6	47,1	49,5
367-1334	10	189,2	167,8	35,1	38,6	42,1	45,6	43,1	45,6	48,1	50,5




Expersus Femoral Stem STD, 135°, cementless (ISO 5832-3 Ti6Al4V) with Head Ø32

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø32/Neck component selected				Stem Offset [mm] when Head Ø32 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1326	2	139,8	118,1	34,6	38,6	42,6	46,6	36,0	38,8	41,6	44,5
367-1327	3	145,0	123,2	35,1	39,1	43,1	47,1	37,2	40,0	42,8	45,7
367-1328	4	150,1	128,6	34,6	38,6	42,6	46,6	38,0	40,8	43,6	46,5
367-1329	5	154,0	132,6	34,6	38,6	42,6	46,6	38,6	41,4	44,2	47,1
367-1330	6	159,9	138,4	34,6	38,6	42,6	46,6	39,0	41,8	44,6	47,5
367-1331	7	164,9	143,5	34,6	38,6	42,6	46,6	39,7	42,5	45,3	48,2
367-1332	8	170,1	148,4	34,6	38,6	42,6	46,6	40,7	43,5	46,3	49,2
367-1333	9	180,0	158,6	34,6	38,6	42,6	46,6	41,8	44,6	47,4	50,3
367-1334	10	189,2	167,8	34,6	38,6	42,6	46,6	42,8	45,6	48,4	51,3




Expersus Femoral Stem STD, 125°, cemented (ISO 5832-12 CoCrMo) with Head Ø28

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø28/Neck component selected				Stem Offset [mm] when Head Ø28 mm/Neck Component selected			
				S (-4)	M (0)	L (+4)	XL (+8)	S (-4)	M (0)	L (+4)	XL (+8)
367-1400	2	138,8	117,6	34,8	38,3	41,8	45,3	41,5	43,5	45,5	47,5
367-1401	3	144,0	122,7	35,8	39,3	42,8	46,3	43,5	45,5	47,5	49,5
367-1402	4	149,1	128,1	34,9	38,4	41,9	45,4	43,9	45,9	47,9	49,9
367-1403	5	153,1	132,1	35,2	38,7	42,2	45,7	44,7	46,7	48,7	50,7
367-1404	6	158,9	137,9	35,1	38,6	42,1	45,6	45,1	47,1	49,1	51,1
367-1405	7	163,9	143,0	35,2	38,7	42,2	45,7	45,9	47,9	49,9	51,9
367-1406	8	169,1	148,2	35,3	38,8	42,3	45,8	47,1	49,1	51,1	53,1
367-1407	9	179,0	158,1	35,4	38,9	42,4	45,9	48,4	50,4	52,4	54,4
367-1408	10	188,2	170,6	35,7	39,2	42,7	46,2	49,6	51,6	53,6	55,6




Expersus Femoral Stem STD, 125°, cemented (ISO 5832-12 CoCrMo) with Head Ø32

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø32/Neck component selected				Stem Offset [mm] when Head Ø32 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1400	2	138,8	117,6	34,3	38,3	42,3	46,3	41,2	43,5	45,8	48,1
367-1401	3	144,0	122,7	35,3	39,3	43,3	47,3	43,2	45,5	47,8	50,1
367-1402	4	149,1	128,1	34,4	38,4	42,4	46,4	43,6	45,9	48,2	50,5
367-1403	5	153,1	132,1	34,7	38,7	42,7	46,7	44,4	46,7	49,0	51,3
367-1404	6	158,9	137,9	34,6	38,6	42,6	46,6	44,8	47,1	49,4	51,7
367-1405	7	163,9	143,0	34,7	38,7	42,7	46,7	45,6	47,9	50,2	52,5
367-1406	8	169,1	148,2	34,8	38,8	42,8	46,8	46,8	49,1	51,4	53,7
367-1407	9	179,0	158,1	34,9	38,9	42,9	46,9	48,1	50,4	52,7	55,0
367-1408	10	188,2	170,6	35,2	39,2	43,2	47,2	49,3	51,6	53,9	56,2



Expersus Femoral Stem STD, 135°, cemented (ISO 5832-12 CoCrMo) with Head Ø28

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø28/Neck component selected				Stem Offset [mm] when Head Ø28 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1317	2	138,8	117,6	35,4	38,9	42,4	45,9	36,3	38,8	41,3	43,7
367-1318	3	144,0	122,7	35,9	39,4	42,9	46,4	37,5	40,0	42,5	44,9
367-1319	4	149,1	128,1	35,4	38,9	42,4	45,9	38,3	40,8	43,3	45,7
367-1320	5	153,0	132,1	35,4	38,9	42,4	45,9	38,9	41,4	43,9	46,3
367-1321	6	158,9	137,9	35,4	38,9	42,4	45,9	39,3	41,8	44,3	46,7
367-1322	7	163,9	143,0	35,4	38,9	42,4	45,9	40,0	42,5	45,0	47,4
367-1323	8	169,1	148,2	35,4	38,9	42,4	45,9	41,0	43,5	46,0	48,4
367-1324	9	179,0	158,1	35,4	38,9	42,4	45,9	42,1	44,6	47,1	49,5
367-1325	10	188,2	167,4	35,4	38,9	42,4	45,9	43,1	45,6	48,1	50,5



Expersus Femoral Stem STD, 135°, cemented (ISO 5832-12 CoCrMo) with Head Ø32

Art.-No.	Size	Length L [mm]	Length L1 [mm]	Neck length [mm] when Head Ø28/Neck component selected				Stem Offset [mm] when Head Ø28 mm/Neck Component selected			
				S (-3,5)	M (0)	L (+3,5)	XL (+7)	S (-3,5)	M (0)	L (+3,5)	XL (+7)
367-1317	2	138,8	117,6	34,9	38,9	42,9	46,9	36,0	38,8	41,6	44,5
367-1318	3	144,0	122,7	35,4	39,4	43,4	47,4	37,2	40,0	42,8	45,7
367-1319	4	149,1	128,1	34,9	38,9	42,9	46,9	38,0	40,8	43,6	46,5
367-1320	5	153,0	132,1	34,9	38,9	42,9	46,9	38,6	41,4	44,2	47,1
367-1321	6	158,9	137,9	34,9	38,9	42,9	46,9	39,0	41,8	44,6	47,5
367-1322	7	163,9	143,0	34,9	38,9	42,9	46,9	39,7	42,5	45,3	48,2
367-1323	8	169,1	148,2	34,9	38,9	42,9	46,9	40,7	43,5	46,3	49,2
367-1324	9	179,0	158,1	34,9	38,9	42,9	46,9	41,8	44,6	47,4	50,3
367-1325	10	188,2	167,4	34,9	38,9	42,9	46,9	42,8	45,6	48,4	51,3

Manufacturer:
OHST Medizintechnik AG
Grünauer Fenn 3
14712 Rathenow
Germany

Fon +49 (0) 3385 5420 0
Fax +49 (0) 3385 5420 99
E-Mail info@ohst.de
URL www.ohst.de

Disclaimer

These informations are intended exclusively for specialists in medical technology. It is not meant as information for medical lay persons.
The explanations about the products contained in the information are of general nature and do not constitute medical advice.

The information was created and compiled by medical experts and technically qualified employees of OHST AG to the best of their knowledge.

No liability or guarantee is assumed for the timeliness, accuracy and completeness of the information provided.
Any liability for material or immaterial damages caused by the use of this information is excluded.

 0482